

## PRODUCT DATASHEET

Model: **WDDR17-C**


Description: **WEIGHT DISTRIBUTION CASSETTE CAP WITH CLEAR TOP FOR 17x17 (43x43cm) DR PANELS**

Revision: April 2018

**Reina**Imaging

### PRODUCT ATTRIBUTES:

- Protects DR Panel during weight-bearing exams.
- Fits over the DR Panel with clear top.
- Features slip resistant grips for stability.
- For non-gridded weight-bearing studies.
- Portable.



### MATERIAL COMPOSITION

PART	MATERIAL	APPROXIMATE %
Shell	Polycarbonate and Polypropylene	98
Overlay	Polycarbonate Film	2

### DIMENSIONS AND PERFORMANCE

Product Dimension: 20.125 x 20.125 x 2.25 in. (51.1 x 51.1 x 5.7 cm)

Weight: 8.4 lbs. (3.81 kg)

Weight Tolerance: 750 lbs. (340 kg)

### FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

### BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

### STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein in this document are accurate and representative of the models being described.

### REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

**Reina Imaging**

**6107 Lou Ave.**

**Crystal Lake, IL 60014 USA**

### AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148Milano Via don Gnocchi 24, Italy

+39 02 48706799

### DEVICE FAMILY CLASSIFICATION

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1

6107 West Lou Avenue | Crystal Lake, IL 60014 | (P) 800.752.4918 | (F) 815.356.8270 | [reinaimaging.com](http://reinaimaging.com)