


Keep pace with supply storage and transport. Keep pace with work flow. Keep pace with patient care.

Pace carts are mobile, easy to maneuver, secure, and adaptable to change. Five cart heights, two widths, three drawer sizes, and a variety of accessories let you target supply storage for specific treatments and departments.





Materials

Shell: **HDPE core with aluminum exterior**
Base: **16-gauge steel with powder-coat finish**
Bumpers: **ABS thermoplastic polymer**
Corner extrusions: **Anodized aluminum**
Work surface: **Antimicrobial ABS thermoplastic polymer**
Drawers: **18-gauge steel with powder-coat finish**
Drawer pulls: **Antimicrobial ABS thermoplastic polymer**
Drawer fronts: **ABS thermoplastic polymer**
Drawer glides: **Zinc plated, steel ball bearing**
Casters: **Heavy duty, 5" full-swivel, ball bearing with steel plate and frame**

Dimensions (nominal)

Exterior: **29.25"w x 21.25"d x 39.25"h**
Weight: **135 lbs**

Options

- 5 case colors and 11 drawer colors
- Keyed, keypad, proximity, and dual access locks
- Back- and side-mounted accessories
- Drawer organization accessories

