


Pace Pediatric Care Cart (SP27PED)

Keep pace with supply storage and transport. Keep pace with work flow. Keep pace with patient care.

Pace carts are mobile, easy to maneuver, secure, and adaptable to change. The Pediatric Care Cart is 29.25" wide and has 9 3"-high, multicolored drawers. A variety of exterior and drawer accessories lets you target supply storage for specific treatments and departments.



Pace Pediatric Care Cart (SP27PED)



Materials

Shell: **HDPE core with aluminum exterior**

Base: **16-gauge steel with powder-coat finish**

Bumpers: **ABS thermoplastic polymer**

Corner extrusions: **Anodized aluminum**

Work surface: **Antimicrobial ABS thermoplastic polymer**

Drawers: **18-gauge steel with powder-coat finish**

Drawer pulls: **Antimicrobial ABS thermoplastic polymer**

Drawer fronts: **ABS thermoplastic polymer**

Drawer glides: **Zinc plated, steel ball bearing**

Casters: **Heavy duty, 5" full-swivel, ball bearing with steel plate and frame**

Dimensions (nominal)

Exterior: **29.25"w x 21.25"d x 39.25"h**

Weight: **135 lbs**

Options

- Keyed, keypad, proximity, and dual access locks
- Back- and side-mounted accessories
- Drawer organization accessories