



## EcoVue® Ultrasound Gel



With a 99% natural formulation free of dyes, parabens and propylene glycol, EcoVue® Ultrasound Gel is gentler on patients and delivers exceptional glide and image quality with less mess and waste. Available in sterile and non-sterile options with a variety of packaging choices, technicians can choose the right amount for the procedure, dramatically reducing waste.

your single source supplier



# Why EcoVue®?

With a 99% natural formulation free of dyes, parabens and propylene glycol, EcoVue® Ultrasound Gel is gentler on patients and delivers exceptional glide and image quality with less mess and waste. Available in sterile and non-sterile options with a variety of packaging choices, technicians can choose the right amount for the procedure, dramatically reducing waste.

Trust EcoVue® to deliver the quality and performance that today's ultrasound procedures demand.



## Single Use Packets

EcoVue® is available in a variety of single-use packets in both sterile and non-sterile options. Our single-use packets are easy to use and specifically designed to help reduce the possibility of cross-contamination associated with reusing bottles. In addition, our single-use packets are also available with a sterile SafeWrap® for an added protection against cross-contamination.



280

**EcoVue® Ultrasound Gel 20g  
Sterile Packet with SafeWrap™**  
48 ea/bx – 6 bx/cs



280NW

**EcoVue® Ultrasound Gel 20g  
Sterile Packet**  
100 ea/bx – 4 bx/cs



281

**EcoVue® Ultrasound Gel 20g  
Non-Sterile Packet**  
100 ea/bx – 4 bx/cs



380

**EcoVue®HV Ultrasound Gel 20g  
Sterile Packet with SafeWrap™**  
48 ea/bx – 6 bx/cs



380NW

**EcoVue®HV Ultrasound Gel 20g  
Sterile Packet**  
100 ea/bx – 4 bx/cs



381

**EcoVue®HV Ultrasound Gel 20g  
Non-Sterile Packet**  
100 ea/bx – 4 bx/cs

- ✓ Choose the right amount of gel for your procedure
- ✓ 99% natural formula is easier on equipment and can reduce wear and tear
- ✓ Provides optimal acoustic velocity at any temperature, even when warmed.
- ✓ Approved by major diagnostic ultrasound equipment manufacturers
- ✓ Cleans up easy with no residue left behind

# Multi-Use Gel

EcoVue® offers a variety of multi-use sizes for diagnostic and therapeutic medical ultrasound procedures and allows technicians can choose the right amount for the procedure.



285

**EcoVue® Ultrasound Gel 250g  
FlexPac® Non-Sterile**  
12 ea/bx – 6 bx/cs



286

**EcoVue® Ultrasound Gel 250mL  
(8.5oz.) Bottle Non-Sterile**  
12 ea/bx – 6 bx/cs



386

**EcoVue®HV Ultrasound Gel 250mL  
(8.5oz.) Bottle Non-Sterile**  
12 ea/bx – 6 bx/cs

## EcoVue®HV High Viscosity Formula

EcoVue's high quality natural formulation is now available in a new high viscosity option, delivering the 99% natural formulation of EcoVue® with a higher viscosity for procedures that demand a thicker gel with even more staying power. Available in a bottle as well as a variety of sterile and non-sterile single-use packets, EcoVue® HV is the latest innovation focused on delivering the quality and performance to meet the changing needs of today's ultrasound procedures.



## Acclimate® Gel Warmer

The patented design allows the device to warm flexible packaged products. The Acclimate® gel warmer is the only warmer designed with single-use packets and standard size flexible packaging in mind. The front “dip” allows for easy access to insert and remove single-use packets without having to reach down into a cavity designed for other types of packaging.

### Models Available:

- **289F** - Acclimate® Gel Warmer - Flexpack® with digital display (1 ea/cs)
- **289B** - Acclimate® Gel Warmer - Bottle with digital display (1 ea/cs)
- **289S** - Acclimate® Gel Warmer - Single-Use Packets with digital display (1 ea/cs)

### More Reasons to Choose Acclimate®

- Wall Mount or Countertop
- Two Stage Thermal Cut-Off Protects from Overheating
- UL Listed
- Available in 120V or 230V
- Two-year Warranty

### Safe Gel Warming for Patient Comfort

The Acclimate® gel warmer provides the convenience of single-button control with a built-in two-stage thermal cut-off to protects from overheating.

Functions as a wall mount or counter top unit, it is UL listed, available in 120V or 230V, and includes a two-year warranty. Acclimate® helps to deliver a better experience for today's ultrasound providers AND their patients.





**Oceania**

PO Box 3225 Loganholme, QLD 4129, Australia  
1300 132 100  
[www.imgsol.com.au](http://www.imgsol.com.au)

**Americas**

222 W. Las Colinas Blvd. Suite 1650E Irving, TX 75039  
800 464 7650  
[www.imgsol.us](http://www.imgsol.us)

**EMEA**

37 Baggot Street Lower, Dublin 2, D02 NV30, Ireland  
+617 3387 0400  
[www.imgsol.eu](http://www.imgsol.eu)

**Asia**

PO Box 3225 Loganholme, QLD 4129, Australia  
+617 3387 0400  
[www.imgsol.asia](http://www.imgsol.asia)

[www.imgsol.com](http://www.imgsol.com) [info@imgsol.com](mailto:info@imgsol.com)

© Imaging Solutions 2022 - Version 0922